



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)


(51) International Patent Classification ⁶ : C07K 14/47, G01N 33/68, 33/53		A1	(11) International Publication Number: WO 98/54219
			(43) International Publication Date: 3 December 1998 (03.12.98)
(21) International Application Number: PCT/US98/10518		(81) Designated States: JP, European patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE).	
(22) International Filing Date: 21 May 1998 (21.05.98)		Published <i>With international search report. Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i>	
(30) Priority Data: 08/865,468 29 May 1997 (29.05.97) US 08/874,566 13 June 1997 (13.06.97) US 08/898,649 22 July 1997 (22.07.97) US			
(71) Applicant: MEDICAL ANALYSIS SYSTEMS INC. (US,US); 542 Flynn Road, Camarillo, CA 94105-93012 (US).			
(72) Inventor: MORIANA, Nihmat; 730 S.W. 99 Avenue, Pembroke Pine, FL 33025 (US).			
(74) Agent: AXFORD, Laurie; Fish & Richardson, Suite 1400, 4225 Executive Square, La Jolla, CA 92037 (US).			

BEST AVAILABLE COPY

(54) Title: COVALENTLY COUPLED TROPONIN COMPLEXES

(57) Abstract

Stable and long-lasting complexes for use in immunoassays for troponin are disclosed. The complexes are covalently coupled complexes of troponin and protein A, and are prepared by reacting troponin with protein A in the presence of a cross-linking agent.

Stability of TnI-153:TnC Complex

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AL	Albania	ES	Spain	LS	Lesotho	SI	Slovenia
AM	Armenia	FI	Finland	LT	Lithuania	SK	Slovakia
AT	Austria	FR	France	LU	Luxembourg	SN	Senegal
AU	Australia	GA	Gabon	LV	Latvia	SZ	Swaziland
AZ	Azerbaijan	GB	United Kingdom	MC	Monaco	TD	Chad
BA	Bosnia and Herzegovina	GE	Georgia	MD	Republic of Moldova	TG	Togo
BB	Barbados	GH	Ghana	MG	Madagascar	TJ	Tajikistan
BE	Belgium	GN	Guinea	MK	The former Yugoslav Republic of Macedonia	TM	Turkmenistan
BF	Burkina Faso	GR	Greece	ML	Mali	TR	Turkey
BG	Bulgaria	HU	Hungary	MN	Mongolia	TT	Trinidad and Tobago
BJ	Benin	IE	Ireland	MR	Mauritania	UA	Ukraine
BR	Brazil	IL	Israel	MI	Malawi	UG	Uganda
BY	Belarus	IS	Iceland	MX	Mexico	US	United States of America
CA	Canada	IT	Italy	NE	Niger	UZ	Uzbekistan
CF	Central African Republic	JP	Japan	NL	Netherlands	VN	Viet Nam
CG	Congo	KE	Kenya	NO	Norway	YU	Yugoslavia
CH	Switzerland	KG	Kyrgyzstan	NZ	New Zealand	ZW	Zimbabwe
CI	Côte d'Ivoire	KP	Democratic People's Republic of Korea	PL	Poland		
CM	Cameroon	KR	Republic of Korea	PT	Portugal		
CN	China	KZ	Kazakhstan	RO	Romania		
CU	Cuba	LC	Saint Lucia	RU	Russian Federation		
CZ	Czech Republic	LI	Liechtenstein	SD	Sudan		
DE	Germany	LK	Sri Lanka	SE	Sweden		
DK	Denmark	LR	Liberia	SG	Singapore		
EE	Estonia						

COVALENTLY COUPLED TROPONIN COMPLEXESField of Invention

5 The present invention relates generally to clinical chemistry. In particular, it relates to stabilized troponin complexes useful in the diagnosis of myocardial infarction or other ischemic events.

Cross-References to Related Applications

10 This application is a continuation-in-part of U.S. Serial No. 08 865,468, filed May 29, 1997 and of U.S. Serial No. 08 874,566, filed on June 13, 1997.

Background of the Invention

15 The determination of the presence or amount of certain constituents or analytes is useful in the diagnosis of disease and physical well-being. Laboratory compositions behave similarly to any constituent found present in human bodily fluids (e.g. blood, blood serum, plasma, spinal fluid or urine). Compositions are used in clinical laboratories. These compositions assist in the determination of whether the clinical
20 instrumentation and procedures used by the laboratory to measure the constituents are accurate. These compositions are also used to calibrate the clinical devices which measure the amount or presence of the constituent in a sample. These compositions will be
25 referred to hereinafter as control compositions or controls.
30

2

In addition, it is important that the analyte or analyte analog present in the control composition behave similarly to the corresponding analyte to be tested for in a patient's bodily fluid --that is, the control composition should
5 mimic the patient sample.

Rapid and simple tests that can be used to accurately diagnose the occurrence of myocardial infarction ("MI" or distinguish other ischemic events
10 such as unstable angina are extremely important. Cardiac troponin I (cTnI) and troponin T have recently become established as the markers of choice in evaluating cardiac distress. See, New England Journal of Medicine Volume 335 No. 19, pages 1342-1349, Antman
15 et al. and pages 1333-1341, Ohman et al.

The Troponin complex is present in both skeletal and cardiac muscles and consists of three subunits, Troponin T ("TnT" the tropomyosin binding subunit, Troponin I ("TnI", the I-actin binding subunit and TnI,
20 which inhibits the actomyosin Mg⁺⁺-ATPase.

The majority of the research into the troponin complex has centered around the regulatory function and structure of the troponin complex in skeletal muscle. The troponin complex assists in muscle contraction.
25 The TnC molecule has four binding domains to bind divalent metal ions. The Ca⁺⁺/Mg⁺⁺ binding sites are in the COOH terminal region and the Ca⁺⁺ binding sites are in the amino terminal region. In studies of skeletal muscle, in the absence of Ca⁺⁺, the amino
30 terminus of TnI binds to the COOH terminus region of

3

TnC and to the globular COOH terminus region of TnT. Thus, research indicates that TnI and TnC are anti-parallel and TnI and TnT are anti-parallel. The presence of calcium ion increases the TnC amino
5 terminus domain's affinity for the inhibitory and COOH regions of TnI. In addition, there is a hydrophobic surface in the N-terminal domain of TnC that represents a Ca^{++} dependent binding site for TnI and TnT. It has been proposed that the Ca^{++} dependent reactions relate
10 to the regulatory mechanism and Ca^{++} independent interactions maintain the structural integrity of the complex. In order to study structure and function of the troponin complex in its regulation of skeletal muscle, cross-linking studies have been accomplished.
15 See, Faran, D. and Reinach, F. Review: The Troponin complex and regulation of muscle contraction. *FASEB Journal*, 9 pp. 763-767, 1995. Covalent binding between TnI and skeletal muscle TnC has been formed between the
20 amino group of lysine in the TnC and lysine groups in TnI using EDC. Miyayoshi et al. 1994, Structure of the troponin complex: implications of photocross-linking of troponin I to troponin C thiol mutants. *J. Biol. Chem.* **269**, 5725-5729. In addition Leszyk et al. (1987) Cross-linking of rabbit skeletal muscle troponin with the
25 photoactive reagent 4-maleimidobenzophenone; identification of residues in troponin I that are close to cystein-98 of troponin C. *Biochemistry* **26**, 7042-7047, reported that the main product of cross-linking between TnC and skeletal muscle TnI comprises segments
30 derived from the N-terminal regulatory domain of TnC

(residues 46 to 78) and the inhibitory region of skeletal TnI (residues 96-116). The Troponin complex is also referred to herein as the ternary complex.

U.S. 08 865,463, filed on May 29, 1997 and also
5 owned by applicant, discloses that it had been discovered that the majority of native cTnI in human serum after an MI is associated with TnC and TnT. The presence of TnI in a complex with other troponin subunits in MI patient serum increases its stability and protects it from further degradation. In addition
10 the troponin complex protects the sites where cardiac-specific antibodies bind. U.S. 08/865,463, filed on May 29, 1997 also discloses methods to isolate the complex from MI patient serum.

15 The cardiac isotype of the myofibrillar contractile protein, Troponin I ("TnI"), is uniquely located in cardiac muscle. TnI is the inhibitory subunit of Troponin, a thin filament regulatory protein complex, which confers calcium sensitivity to the
20 cardiac and striated muscle. Troponin I exists in three isoforms: two skeletal TnI (fast and slow) isoforms (Molecular Weight = 19,800 daltons) and a cardiac TnI ("cTnI") isoform with an additional 31 residues (human TnI) on the N-terminus resulting in a molecular weight
25 of 23,000 daltons.

Cardiac TnI is found in human serum rapidly (within approximately 4 to 6 hours) following an MI. It reaches a peak level after approximately 18-24 hours and remains at elevated levels in the blood stream for
30 up to 6 to 7 days. Thus, immunoassays which can test

for human cTnI are valuable to the medical community and to the public.

It is desirable to use an immunologically reactive human cTnI isoform comparable to that found in MI patient serum. We found that MI patient serum contains TnI fragment(s) which is the result of the C-terminal processing of cTnI molecule. The high sequence homology found in the C-terminal region between cardiac TnI and skeletal muscle TnI (Larue et al. 1992 Molec. Immunology 39, 271-278, Vallins et al. 1990 FEBS Lett. 270, 57-61, Leszky et al. 1988 Biochemistry 27, 2821-2827) produce TnI antibodies directed against this region having non-cardiac specificity (Larue et al. 1992). Our data and Larue et al. 1992 suggest that most of the known cTnI specific antibodies have their epitopes located approximately in the first 75% of the TnI molecule. Therefore, this portion of the TnI molecule should function as a MI specific cTnI isoform in most immunoassay systems.

Currently, cTnI immunoassays are commercially available from Dade International, Behring Diagnostics, and Sanofi Pasteur Diagnostics. The Dade product is the Stratus® Cardiac Troponin-I assay.

Native intact human cTnI is difficult to obtain because of the scarcity of human heart and native intact human cTnI is highly subject to proteolytic degradation during purification. Recombinant cardiac TnI ("r-TnI"), unlike the native human cTnI, can be produced and purified in acceptable quantities. As expressed by Dade, the primary structure of r-TnI

6

contains 226 amino acids (SEQ ID NO: 1); 209 of them represent the TnI sequence (SEQ ID NO: 2). In addition to the primary sequence of cTnI (SEQ ID NO: 2), r-TnI, as expressed by Dade International, has a leading sequence of 8 amino acids (MASMTLWM) on the N-terminal, and a tail sequence of 9 amino acids (PMVHHHHHH) on the C-terminal (SEQ ID NO: 1). The primary structure of the r-TnI molecule has methionine residues at positions -7, -4, 0, 153, 154, 200 and 211 (SEQ ID NO: 1). See also Figure 1.

Full length cardiac troponin I is known to have the following sequence:

```
ADGSSDAAREPRPAPAFIRRRSSNYRAYATEPHAKKKSKISASRKLQLKT
LLIQIAHQBELEREAEERRGKNGRALSTROQPLELTGLGFAELQDLORQLH
RNVQWVDEEPYDIEAKVTKNITEIADLTQKIFDLRGFWFPPILRPVRISA
DAMMIALIGAPAHESLCLPANHLYNNKEDTEKENPEWGIWPKNTIALSGME
DPYFNFES
```

SEQ ID NO: 2 Armour, M.L. et al., 1993 Cloning and Expression in *Escheria coli* of the cDNA Encoding Human Cardiac Troponin I. Gene, 131:1-197-191.

U.S. Serial No.08/564,526, also owned by applicants, and incorporated herein by reference, discloses the use of a human cTnI fragment generated from human r-TnI by chemical cleavage. The cleavage of r-TnI by cyanogen bromide (CNBr) results in a major polypeptide of 153 amino acids, hereinafter referred to as the "CNBr-cTnI isoform" (SEQ ID NO: 3). See Figure 2. The CNBr-cTnI isoform represents 73% of the primary structure of human cTnI and is immunologically more reactive than r-TnI. The purified CNBr-cTnI isoform

7

has an average of 3-4 times more reactivity than r-TnI and lower non-specific binding, as measured by radial partition immunoassay. The molecular size of the CNBr-cTnI isoform is comparable in molecular weight to the major degradation product of native cardiac TnI in MI patient serum.

It is desirable to use an immunologically reactive human cTnI isoform comparable to that detected in MI patient serum. The availability of r-TnI can facilitate the production of cardiac cTnI isoforms. Moreover, since most of the known human cardiac specific TnI antibodies have their epitopes located approximately in the first 75% of the TnI molecule, that portion of the TnI molecule will function as a cTnI isoform in most immunoassays.

The CNBr-cTnI isoform can be used as calibrators or controls in various cTnI immunoassays.

U.S. Pat. 5,453,406, filed May 22, 1997, discloses the use of cardiac troponin I fragments of the general sequence X-A-B-Y wherein X comprises any of amino acids 1-27 of full length cardiac troponin I, A comprises residues 28-69 of full length cardiac troponin I, B comprises amino acid residues 70-90 of full length cardiac troponin I, and Y comprises any sequential amino acid sequence of amino acid residues 91-170 of full length cardiac troponin I. These sequences also have increased immunoreactivity and stability over prior art compounds.

8

Troponin T (TnT) with a molecular weight of 39,000 Kd is part of the troponin-tropomyosin complex of the thin filament that is part of the muscle contractile apparatus and that contains actin and tropomyosin regulatory elements. Skeletal muscle studies of TnT have found that TnT is structurally asymmetric. Its terminal globular COOH terminal domain (TnT-3) mediates its interaction with TnI and TnC. TnT-1 at the amino terminus domain interacts with tropomyosin. See, Farah, C. and Reinach, F. (1995) Review: The Troponin complex and regulation of muscle contraction. *FASEB Journal* 9, 755-767. It is reported that skeletal TnT is cleaved into the skeletal TnT-1 and TnT-2TnI-TnC fragments by mild proteolysis. Schaertl, S. et al. (1995) Separation and Characterization of the Two Functional Regions of Troponin Involved in Muscle Thin Filament Regulation. *Biochemistry* 34, 490, 15991-15994. TnT serves as a link between the tropomyosin backbone and the TnI and TnC troponin complex. TnT has isotypes in cardiac and fast and slow skeletal muscles. It appears in serum about 3 hours after the onset of chest pain and remains elevated for at least 10 days following MI. Despite its lack of complete cardiac specificity it can be useful because of its rapid appearance into the bloodstream. Troponin T can be obtained as described in *J. Biochem.* 72: pages 723-735 (1972) or *J. Biol. Chem.* 249: 4742-4748, or purchased commercially. TnT gene promoter and derivatives thereof are disclosed in U.S. 5,266,433. TnT isoforms of skeletal muscle show variation in a given species in about a 30 amino acid

9

region of the amino terminus and about a 14 amino acid region of the carboxy terminus. Pan, B.S. and Potter, J.D. (1990): Two Genetically Expressed Troponin T Fragments Representing α and β Isoforms Exhibit Functional Differences. *Journal of Biological Chemistry* **267** (82): 23052-23056.

In vitro stabilized solutions for cardiac markers have been disclosed. U.S. 5,533,200 and Bodor et al., (1992) Development of Monoclonal Antibodies for an Assay of Cardiac Troponin-I and Preliminary Results in Suspected Cases of Myocardial Infarction, *Clinical Chemistry* **38**, (11): 2203-2214 at 2204 disclose stabilized troponin T and/or troponin I using troponin C and calcium ion. U.S. 5,533,200 discloses that serum may be added. U.S. 08,874,566, filed June 13, 1997, discloses improvements in stabilizing the troponin T or troponin I troponin T complex and discloses solutions useful as calibrators or controls for diagnostic assay, including troponin. U.S. Serial No. 08/544,114 and U.S. Serial No. 08/545,465, filed May 27, 1997 also disclose the effect of TnC upon the immunological and biological activity and non-specific binding of the CNBr-cTnI isoform and other fragments. U.S. Serial No. 08/564,526 discloses the activity of the complex formed by the CNBr - cTnI isoform, TnC and TnT as useful in immunoassays.

The calibrators and controls in Behring's OPUS® assay are a lyophilized preparation of human cardiac troponin I in processed bovine calf serum with stabilizers. The reconstituted products are stable for

10

seven days when stored at 2 to 8C. The calibrators and controls in Sanofi Pasteur's troponin I assay are a lyophilized preparation in a buffered human serum matrix. The reconstituted calibrators must be used within fifteen minutes after complete reconstitution, but may be aliquoted and stored frozen at -20C for up to about six months. The calibrators and controls in the Dade troponin I assay are provided frozen. When thawed the product is stable for thirty days when stored at 2-8C.

Improved methods and compositions of stabilized troponin I and/or troponin T are still necessary because under certain conditions the complex can dissociate (e.g. removal of calcium, presence of detergents such as SDS). Moreover, the closer the analyte used in the composition is to the actual circulating isoform of the troponin complex, the better the composition will perform as a primary reference material --that is a calibrator on which other calibrators are based.

Summary of the Invention

This invention relates to stabilized compositions of: troponin I and troponin C complex (TnI:TnC); troponin T and troponin C complex (TnT:TnC); troponin I and troponin T complexes (TnI:TnT) and troponin I, troponin T, and troponin C complexes (TnI:TnC:TnT) for use in immunoassay of cardiac troponins.

Troponin I and/or troponin T can be covalently complexed to troponin C to provide a composition that

has enhanced stability and/or immunoreactivity over prior art complexes and analytes.

In a TnI:TnC complex, the TnI is covalently coupled to TnC. The troponin I may be native or recombinant or may be fragmented or full-length.

In a TnT:TnC complex, the TnT is covalently coupled to TnC. The TnT may be native or recombinant or may be fragmented or full-length.

In a TnT:TnI complex the TnI and TnT are covalently coupled.

In a TnT:TnC:TnI complex either or both of the TnI or TnT is covalently coupled to TnC.

The complexes are useful as calibrators or controls for methods that assay for TnI, TnC, and/or TnT or for use as primary reference materials.

The control composition should contain a buffer or serum base matrix and may contain such metal ions as calcium and magnesium ions. The control composition should be provided in liquid form.

The covalent coupling agents include those coupling agents that provide substantially "native-length" covalent cross-linking between troponin T and/or troponin I with troponin C. The term "native-length" cross-linking as used herein means a covalent bond formed between either troponin I or T with troponin C that provides a covalently coupled complex that has substantially the same immunological activity as non-covalently coupled complex. Generally the length of the covalent bond should approximate the length between troponin I and/or T with troponin C in

12

native complex or between troponin I and troponin T in native complex. But, because protein complexes are rarely rigid in structure, it is to be understood that there is variability in the structure of the complex.

5 The stability and immunological activity of the covalently coupled complex are what is important.

Brief Description of the Figures

Figure 1 shows the amino acid sequence of cardiac troponin as expressed by Dade International (SEQ ID NO. 1) and of cardiac troponin I (SEQ ID NO. 2).

5 Figure 2 shows the amino acid sequence of the CNBr cleavage product of recombinant troponin I (human cTnI isoform or TnI-153).

Figure 3 shows the stability and immunoreactivity at 45 C of cross-linked rTnI-153:TnC complex at two
10 dilutions of complex. The samples were evaluated using a Stratus II analyzer and Stratus Troponin Fluorometric Assay Kit both available from Dade International.

Figure 4 shows the stability at different temperatures of diluted cross-linked rTnI-153:TnC
15 complex. The samples were evaluated using a Stratus II analyzer and Stratus Troponin Fluorometric Assay Kit both available from Dade International.

Figure 5 shows the stability of the cross-linked rTnI-153:TnC complex using three different lots of TnC.

20 Figure 6 shows the stability of cross-linked rTnI:TnC complex at 45 C.

Figure 7 shows polyacrylamide gels of cross-linked and non-crosslinked complexes of rTnI-153:TnC.

Detailed Description of the Invention

In a TnI:TnC complex, the TnI is covalently coupled to TnC. The cardiac troponin I may be native or recombinant and may be fragmented or full-length. While some uncomplexed troponins can be found in human serum after a myocardial event, most cardiac specific troponin is found as complex. It has been found that the TnI in the complex is degraded by proteolytic cleavage at the C-terminal end to provide an 18,000 Kd fragment and a 14,000 Kd fragment. Generally the 14,000 Kd fragment is cleaved from the 18,000 Kd fragment. After the cleavage to the 18,000 fragment, an N-terminus proteolytic cleavage occurs at the carboxyl side of Arg 26, thus eliminating the first 26 amino acids of the N-terminus. The 31 amino acid sequence at the N-terminus had been proposed as the best position to direct antibodies against, however, this recent finding suggests that these antibodies would only recognize a fraction of serum TnI.

Thus, it is preferred that the complexes of the present invention that contain cardiac troponin I contain at least a fragment of troponin I generated from the 14,000 Kd fragment. It is inherent that antibodies for use in immunoassay be generated against that portion of troponin I or troponin complex that includes the 14,000 Kd troponin I sequence. Of course, the antibody must react immunologically (e.g. have an epitopic site on the fragment) and specifically (e.g. it should not substantially cross react with skeletal

15

muscle troponin I) with the cardiac troponin I or complex used in the calibrator or control.

Thus, the use of a human cTnI fragment generated from human cardiac r-TnI by chemical cleavage is a preferred fragment for the TnI:TnC complex because it is closest to the majority of the native form. The cleavage of r-TnI by cyanogen bromide (CNBr) results in a major polypeptide of 183 amino acids, hereinafter referred to as the "CNBr-cTnI isoform" (SEQ ID NO: 3). The CNBr-cTnI isoform represents 73% of the primary structure of human cTnI and is immunologically more reactive than r-TnI as determined using radial partition immunoassay. The purified CNBr-cTnI isoform has an average of 3-4 times more reactivity than r-TnI and lower non-specific binding, as measured by radial partition immunoassay, available from Cade International Inc. As demonstrated in Figure 1 the molecular size of the CNBr-cTnI isoform is comparable in molecular weight to the major degradation product of native cardiac TnI in MI patient serum and retained the epitopes for the antibodies used in the Stratus[®] II TnI Immunoassay System. (See Vallins et al. (1990) *FEBS Lett.* 270, 57-61.)

Generally described, the first step in cyanogen bromide cleavage is to carboxymethylate the cysteine residues of r-TnI (there are two in the TnI sequence) (SEQ ID NO: 1) at positions 79 and 96 in order to prevent dimerization by inter or intra molecular disulfide bridges. The carboxymethylation of the cysteine residues is not a pre-requisite for the

16

generation of the 153 amino acid isoform. (SEQ ID NO: 3.) Rather, the carboxymethylation facilitates the process by minimizing the complications during or after CNBr digestion.

5 CNBr treatment is carried out on the carboxymethylated r-TnI. Unlike other possible cleavage reactions (e.g. enzymatic), the CNBr treatment removes the tail sequence, the leading sequence, and part of the TnI C-terminal region without affecting the
10 primary sequence of the immunogenic sites.

Other preferred peptides are disclosed in U.S. 08/368,468, filed May 29, 1997 and include cardiac troponin I fragments of the general sequence X-A-B-Y wherein X comprises any of amino acids 1-27 of full
15 length cardiac troponin I, A comprises residues 28-69 of full length cardiac troponin I, B comprises amino acid residues 70-90 of full length cardiac troponin I, and Y comprises any sequential amino acid sequence of amino acid residues 91-118 of full length cardiac
20 troponin I.

Preferred residues for X include residues 1-27, 2-27, 3-27, 4-27, 5-27, 6-27, 7-27, 8-27, 9-27, 10-27, 15-27, 20-27, 21-27, 22-27, 23-27, 24-27, 25-27, 26-27, and 27 of SEQ ID NO: 2. More preferably, X is amino
25 acid 27 of SEQ ID NO: 2.

A comprises amino acid residues 28-69 of SEQ ID NO: 2. B comprises amino acid residues 70-90 of SEQ ID NO: 2. Preferred residues of Y include residues 91-92, 91-93, 91-94, 91-95, 91-96, 91-97, 91-98, 91-99,
30 91-100, 91-105, 91-110, 91-115, 91-116, 91-117, 91-118,

17

91-119, 91-120, 91-121, 91-122, 91-123, 91-124, 91-125, 91-126, 91-127, 91-128, 91-129, 91-130, 91-131, 91-132, 91-133, 91-134, 91-135, 91-136, 91-137, 91-138, 91-139, 91-140, 91-141, 91-142, 91-143, 91-144, 91-145, 91-146, 91-147, 91-148, 91-149, 91-150, 91-151, 91-152, 91-153, 91-154, 91-155, 91-160, 91-165, 91-170 of SEQ ID NO: 2. More preferably, Y can be any of 91-95, 91-100, 91-105, 91-110, 91-115, 91-120, 91-130, 91-140, 91-145, 91-150, 91-153, 91-155, 91-160, 91-165, 91-170 of SEQ ID NO: 2.

The lower molecular weight 14,000 TnI fragment, isolated from a pool of patient serum has been sequenced for N-terminal identification. The N-terminal sequence of the TnI 14,000 fragment starts at position 17 Ala in human cardiac TnI sequence. The 14,000 fragment is approximately 111 amino acids long, ending in the region from about amino acid 120 to about amino acid 131 in intact cTnI. The N-terminal sequence of the 14,000 fragment starts at a very near the N-terminus of the intact human cTnI. The 15,000 fragment is approximately 140 amino acids long, ending in the region from about amino acid number 135 to about 145 in intact cTnI. Thus, one preferred group of fragments has X as 25-27, 26-27 or 27 of SEQ ID NO: 2 and Y as 91 to any of 135-145 of SEQ ID NO: 2.

The fragments also may be those cardiac troponin I protein fragments containing the sequence AYATEPHAKKKSKISASRKLQLKTLLLQIAKQEL (SEQ ID NO: 4) or RAYATEPHAKKKSKISASRKLQLKTLLLQIAKQEL (SEQ ID NO: 5). The fragments may be recombinant sequences such as

18

MADGSSDAAREPRPAPAPIRRRSSNYRAYATEPHAKKKSKISASRKLQLKTLLLQ
IAKQELEREAEEREGEKGRLSTROQPLELAGLGFAELQDLQRQLHARVDKVDDEE
RYDIEAKVTKNITEIADLTQKIFDLGKFKRPTLRRVRISADAMMQALLGARAKE
SLDLRAHLKQVKNEDTEKENREVGDWKKNIDALSGMEGRKXKKKFEES (SEQ
ID NO: 6);

ADGSSDAAREPRPAPAPIRRRSSNYRAYATEPHAKKKSKISASRKLQLKTLLLQI
AKQELEREAEEREGEKGRLSTROQ (SEQ ID NO: 7); or similar
human or bovine fragments. It should be understood
that while human troponin is preferred, other species
may be substituted. Generally, these other species
lack the appropriate methionine residue in the full
length primary structure. However, insertion of
methionine into the primary structure at positions that
upon cleavage by CNBr would provide an appropriate
fragment would allow the use of alternative species.

Troponin C is commercially available from a number
of sources and the source or species is not critical.
Generally, rabbit troponin C is used because of the
lower cost, but human and other species can also be
used. Recombinant troponin C can also be used.
The molecular weight of troponin C is about 17,500 kD.
As stated earlier, TnC has a $\text{Ca}^{++}/\text{Mg}^{++}$ binding domain
in the COOH terminal region and a Ca^{++} binding domain
in the amino terminus region and is thought to be
"dumb-bell" shaped connected by a long central helix.

The covalent coupling agents useful in the present
invention include those coupling agents that provide
substantially "native-length" covalent cross-linking
between troponin T and/or troponin I with troponin C.
The preferred coupling agents form covalent bonds

between the TnI or TnT and TnC that approximate the distance between the proteins as they exist in the complex. The cross-linkers may be those "zero-length" crosslinkers such as 1-ethyl-3-(3-(dimethylamino) propyl) carbodiimide (EDC) and other water soluble carbodiimides providing an amide bond between an activated carboxyl group on one protein and coupling to an amino group on the other protein. See, Chemistry of Protein Conjugation and Cross-Linking, CRC Press, Inc.

Other coupling agents such as oxalic, malonic and succinic acid derivatives may also be used as may many of those commercially available cross-linkers such as those available from Pierce Chemical Co. The actual coupling agent used is less important than how the covalently coupled complex performs immunologically or with regard to stability when compared to non-covalently coupled complex, native complex or commercially available products.

The complex is prepared by reacting the coupling agent with troponin C and troponin I in a manner specific for the coupling agent. The coupling reaction of the TnC and TnI should take place in the presence of Ca^{++} or Mg^{++} or other divalent metal ions. The conditions will vary depending on the coupling agent utilized, but those conditions can readily be determined by those skilled in the art with reference to known and published methods. The choice of buffer is not critical, although the buffer should have a molarity from about 10 to 200 mM. The concentration of calcium and magnesium ion is not critical, but

20

preferably should be from about 20 μ M to about 20 mM.
If magnesium ion is used alone it should be in higher
concentrations than that used for calcium. Typical
amounts of calcium and/or magnesium are from about 2-5
5 mM. The buffer may include a salt such as sodium
chloride or potassium chloride at from about 50 mM to
500 mM. The amount of troponin C is not critical, but
should be in molar amounts equal to or greater than the
amount of troponin I that will be utilized. Generally,
10 the amount of troponin C may be from 0.02 mg/mL to 5
mg/mL. The amount of coupling agent may need to be
adjusted depending on the amount of troponin C.

Clinically significant concentrations of diluted
final coupled product are evaluated immunologically in
15 a troponin I assay (such as those available from Dade
International Inc. and may be compared to native
complex or non-covalently coupled complex or
commercially available products evaluated in the same
manner. Generally the clinical range of interest of
20 troponin I or T is about 1.11 ng/mL or less than 1
 μ g/mL and typically 0.1 ng/mL to 200 ng/mL. In
addition, stability may be determined by comparing the
immunological activity over time at one or more
temperatures in comparison to native complex or non-
25 covalently coupled complex or commercially available
products evaluated in the same manner or by running
electrophoretic gels of the complexes to evaluate the
stability.

Native complex may be obtained in serum samples.
30 In addition, native complex can be isolated. Methods

are known to those skilled in the art. For instance, a preferred method of isolating the Tn complex comprises incubating the sample to be tested with a substrate coated with antibodies to the subunits of the Tn complex. Many antibodies are useful, and can be selected by one of skill in the art. Examples of such antibodies include anti-TnI, anti-TnC, and anti-TnT antibodies. A preferred substrate is beads, and more preferably the substrate comprises latex beads. The bound complex is eluted under conditions that do not affect association of Tn subunits, e.g. using urea. These conditions can be determined by one of skill in the art. A preferred buffer system comprises urea and lacks SDS.

A preferred coupling agent is EDC or other water soluble carbodiimides. In a buffered solution containing salt and calcium ions, troponin C is reacted with EDC. The concentration of EDC is not critical and may be from 1-1 mM. Compounds such as N-hydroxysuccinimide (NHS) or SNHS at about 1-1 mM may be added to enhance the activation of troponin C by water soluble carbodiimides. The NHS is typically added prior to the addition of EDC and incubates with the protein for at least about 5 minutes, although 15 minutes is typical. Then the EDC is added and the mixture is incubated, typically at room temperature, for about 15 minutes, but typically for thirty minutes or more. The reaction proceeds best at mildly acidic pH values (e.g. about 6), but the pH may range from 5-9. After a sufficient reaction time, a reducing agent

such as mercaptoethanol is added to stop the reaction. Other reducing agents may be substituted for mercaptoethanol.

Next, cardiac troponin I, preferably the CNBr-TnI isoform, is added to provide a molar ratio of about 1:1 with troponin C, although the amount of TnI may be less. TnI is generally in a buffer that contains salt or urea in sufficient amounts to solubilize the TnI. A typical buffer is 100 mM sodium phosphate, 10 mM Tris and 3 M urea at pH 8 (PTU buffer). The choice of the buffer is not critical to the invention. The buffer, however, must maintain the solubility of the cTnI. The reactants incubate, generally at room temperature, for about an hour and typically about two hours. The formed complex may be buffer exchanged into another buffer. The activity of the cross-linked complex is measured using a Stratus Troponin immunoassay and compared to native complex. The presence of the covalently cross-linked complex may be confirmed using polyacrylamide gel electrophoresis. The complex will not dissociate under reducing conditions or in the presence of EDTA or other metal complexing agent.

When using the non-specific coupling agents such as EDC, calcium or magnesium ion should be present in either the TnC or TnI solution. The calcium and/or magnesium causes complex formation to occur. Thus, the complex forms as it would in native complex so the spatial orientation of the proteins will be comparable to native complex or non-covalently coupled complex. Thus, the cross-linking occurs more selectively and the

23

bond formed will "freeze" the complex in the correct conformation. The covalently coupled complex is insensitive to EDTA or other chelating agents (thus, independent of calcium or magnesium), does not
5 dissociate in the presence of SDS, and is more resistant to higher temperatures than the non-cross-linked complex.

In addition, generally the TnC, rather than the TnI or TnT, should be subjected to EDC or other water
10 soluble carbodiimide treatments. For instance, if TnI is activated by EDC without protection of the epitopic site, there is the possibility that the epitopic site would be affected by the non-specific coupling agent -- that is the TnI or fragment would be less
15 immunologically active.

The TnT:TnC complex and TnI:TnC:TnI complex are prepared in a similar fashion to the TnI:TnC. TnT fragments, particularly the carboxy fragments such as TnT-1, are protected. As described for troponin I, the
20 troponin T should be added in molar amounts equal to or less than the amount of troponin C. Troponin C may be from less than 0.02 mg/mL to 5 mg/mL. The ranges for TnI and TnT are similar. The TnT:TnC complex is evaluated as described for troponin I complexes, except
25 using an assay for troponin T. Troponin T assays are available from Boehringer Mannheim.

A TnT:TnI complex may also be prepared, however, as with the other complexes of this invention, the
immunological (epitopic) sites must be substantially
30 retained. Since both TnT and TnI are measured

immunologically, it is preferred that the epitopic region be retained for both proteins. Thus, less non-specific coupling agents may be useful than agents such as EDC or oxalic acid derivatives.

5 It has been described that skeletal muscle TnI complexes with TnT at about amino acids 40-80 of TnI. See Potter, J.D. et al. (1995) A Direct Regulatory Role of Troponin T and a Dual Role for Troponin C in the Ca++ Regulation of Muscle Control. *Journal of Biological Chemistry* **270** (6) 2557-2562. Thus, it is proposed that the homologous site on cTnI would be useful for coupling. Preferably, TnI and TnT are added in approximately equimolar ratios. The amounts may be from 0.02 mg/mL to 5 mg/mL. It is important
10 that the immunological site not be altered. It should be understood that the epitopic site will vary depending on the antibody used in the assay. The generation of antibodies is described in the art, for example, by Bior et al., 1992 Development of Monoclonal Antibodies for an Assay of Cardiac Troponin-I and Preliminary Results in Suspected Cases of Myocardial Infarction, *Clinical Chemistry* **38** (11) 2203-2214. In addition, several methods that assay for troponin I or T are commercially available.
15

25 The present invention uses a defined base material useful as a control matrix for containing and maintaining the troponin complexes, the base material comprising an aqueous solution of a buffer to maintain the pH at 5-8, anti-microbial

agents and may contain other stabilizers including calcium ions or proteins.

The defined base material may be utilized to prepare stock solutions and controls including troponin I or troponin I fragments, troponin T or troponin T fragments complexed with Troponin C or troponin I complexed with troponin T.

The resulting solution can be stored as a liquid or frozen or also can be lyophilized if appropriate fillers are included.

It is preferred to prepare a stock solution at higher concentrations of troponin I and/or troponin T than those that will be utilized in the final assay control. The stock can be stored frozen or lyophilized and thawed or reconstituted when necessary to prepare the appropriate dilutions of controls or calibration standards. Thus, the complexes prepared as described above may be further exchanged into an appropriate buffer after covalent coupling.

The aqueous solution used to prepare the stock and controls may include a buffer and the buffer may generally be any of the buffers that function in the pH range of 5 to 8. Of these buffers, the buffers that are preferred function are in the pH range of 6-8. The concentration of buffer is between 10 mM to 200 mM. It is preferred to keep the buffer concentration lower - in the range of 20-100 mM. Preferably, the buffers used for the calibrators or controls contain bovine serum albumin (BSA). In certain embodiments the buffer

contains, BSA, sugars, salt and an antibacterial agent. Examples of useful buffers include HEPES, MES or TRIS buffers. A preferred buffer comprises MES buffer containing 6.5% BSA, at pH 6-7. Other preferred
5 buffers contain a reducing agent, stabilizing protein, chelating agent and a salt as described in the copending application U.S.S.N. 08/400,153, incorporated herein by reference.

Alternatively, instead of a buffer, the
10 fragments are spiked into the serum, e.g., human or bovine, or into diluted serum, e.g., serum diluted 1:1 with MES buffer containing BSA.

Anti-microbial and anti-fungal agents may be added to prevent growth and may include those
15 commonly found in the prior art at the concentrations found in the prior art such as pentamycin, clostrimazole, sodium azide, nystatin, amphotericin, Vioform and or Fungilin 314.

In addition, stabilizing proteins such as
20 albumin, gelatin, ovalbumin, or casein may be included. The concentration of stabilizing protein may be from 0 to 15% and preferably from 7 to 12%. Preferably the stabilizing protein is albumin and preferably the albumin is substantially protease
25 free.

It is preferred that the solution have low protease activity, thus protease inhibitors such as aprotinin and "Protease Inhibitor" (Sigma) are
30 effective. However, the use of the recombinant fragments as described herein are not as sensitive

27

to protease activity as is the full length protein. The inhibitors may be added and may be used at the manufacturer's recommended concentration.

5 Examples of other protease inhibitors include benzamidine, (1S, 3R)-3-Amino-2-hydroxy-5-methylhexanoyl]-Val-Val-Asp (Amastatin-Sigma), [2S,3R]-3-Amino-2-hydroxy-4-[4-nitrophenyl]-butanoyl-L-leucine, Antipain, [2S,3R]-3-Amino-2-hydroxy-5-methylhexanoyl]-Val-Val-Asp
10 (Epiamastatin-Sigma), ([2R,3R]-3-Amino-2-hydroxy-4-phenylbutanoyl)-L-leucine (EpiBESTATIN-Sigma), Foroxymithine, Acetyl-Leu-Leu-Arg-al (Leupeptin-Sigma), 4-Amino-3-hydroxy-6-methyl-heptanoic acid,
15 4-Amino-1-hydroxy-6-methylheptanoic acid, N- α -Bammapyrrolysyl-L-phenylalanine (PMSF). It is most preferred that the means to provide a substantially protease free solution is to use
20 substantially protease free proteins such as albumin which is substantially protease free.

25 Serum may be included if desired. Again, the use of fragments that are similar to the 14,000 Kd fragments substantially eliminates the concerns of proteases.

30 Controls prepared by this method may be lyophilized by adding those bulking agents that are known in the art, but the controls may also be liquid. In addition, the liquid controls may be frozen to further increase shelf-life.

Currently used assays for detecting cTnI and cTnT in MI patient serum utilize a sandwich assay. However, the complexes of the present invention can also be used to design competitive-type assays for the detection of cTnI or cTnT in serum. In such an assay, a subsaturating amount of antibody to cTnI or cTnT is bound to a solid phase, e.g., a microtiter plate or latex beads. The complexes of the present invention are labeled, e.g., with alkaline phosphatase, or horseradish peroxidase. A constant amount of the labeled complex is mixed with the sample of MI patient serum containing an unknown amount of cTnI and/or cTnT. The test sample is then allowed to bind to the subsaturating amount of cTnI and/or cTnT antibody bound to a solid phase. The cTnI and/or cTnT in the sample will compete with the labeled complex for binding with the antibody-coated solid phase. Unbound proteins are removed by washing and the amount of labeled complex bound to the solid phase is measured. The amount of labeled complex bound to the antibody on the solid phase indicates the amount of cTnI or cTnT present in the serum. If the serum contains a high concentration of cTnI or cTnT, it will compete effectively with the labeled complex and little or none of the labeled complex will bind the antibody-coated solid phase.

Changes in some amino acids of the fragments of TnI, TnT, and TnC might not affect their performance except those occurring at the epitope(s) where the specific assay antibodies bind and those amino acids of the binding domains for TnC. The affects can readily

be determined by gel electrophoresis and immunological analysis of prepared complex.

Example 1

5 Preparation of a defined base material for cardiac markers.

Antioxidants such as 200 milligrams of glutathione, 200 milligrams of glucose, 50 mgs of ascorbic acid, and 1.1 milliliter of phenol, about 10 2.7 grams of L-Lactate, about 225 milligrams of calcium chloride or other calcium salt to provide 1-3 mM calcium ion, anti-microbial and anti-fungal agents such as about 20 milligrams of chlortrimazole, 35 milligrams of gentamicin, and 15 about 1 milliliter of Proclin 300, about 25 grams of protease free BSA, and about 1 gram of gelatin are combined in an aqueous 50 mM TRIS buffered solution at about pH 7.3 also containing a salt such as sodium chloride about 10 grams to provide about one liter of base material.

It is best to add the gelatin in solution by dissolving the gelatin by adding 1 gram of gelatin to 100 milliliters of water and gently heating to dissolve the gelatin. Then the gelatin containing 25 solution is added to the base material.

The resulting solution is filtered using filters sufficient to remove any bacteria such as a .22 micron filter. A low protein binding filter is preferred.

Example 2

Cross-Linking of Troponin I and Troponin TPreparation of troponin I stock solution and controls.

5 The cross-linking was accomplished using the EDC method described by Greaser and Gergely (1971). To 20 mM Mes buffer at pH 6 containing 0.1 M KCl, 1.2 mM CaCl₂, freshly prepared NHS at 5.7 mM final concentration was added followed by addition of TnC to provide 0.2 mg/mL. The mixture was incubated at room temperature for 15 minutes. The TnC was activated by addition of 5.7 mM EDC (final concentration) followed by incubation at room temperature for 30 minutes. The activation step was terminated by addition of β mercaptoethanol (20 mM final concentration). The CNBr-rTnI isoform TnI 130 in PTU buffer was added at a final concentration of 0.175 mg/mL. The mixture was incubated for two hours at room temperature. The complex is buffer exchanged. The activity of the crossed linked complex was measured using Stratus TnI immunoassay. The presence of a covalent cross-linked complex was confirmed on polyacrylamide gel electrophoresis. The resulting complex is less sensitive to environment, including temperature than native complex. The resulting stock solution is sterile filtered using low protein binding filters of 0.22 microns or less. The stock solution may be stored frozen at -20° C for longer than two years. The

stock solution may be stored at 2-8° C for more than one week.

The stock solution is used to prepare diluted solutions of troponin I in the clinical range of interest (about 0 to 200 ng/mL) by diluting the stock solution in the defined base material prepared in Example 1 or other base material.

Example 3

Stability and Immunological Activity of the Covalently Coupled Complex.

Various aliquots representing different concentrations of troponin I complex prepared in Example 2 of the diluted stock solution, were stored at various temperatures. The aliquots were analysed on a Stratus II Immunoassay analyzer for troponin I concentration. The change in concentration was evaluated over time. Figure 3 shows the stability of two concentrations of diluted stock (15 ng/mL and 45 ng/mL of troponin I) of the complex stored at 45 C. The dilutions are stable for greater than three weeks at 45 C. Figure 4 shows a 28 ng/mL solution was evaluated for stability at 4C, 25C, 37C, and 45 C and was stable for at least six weeks at all temperatures. Figure 5 shows the stability of the complex when prepared with three different lot numbers of commercially available TnC. All lots were stable for over 10 days at 45 C. A full length rTnI-TnC

complex was also prepared. The stability of the cross-linked complex at 45 C is shown in Figure 6.

Example 4

5 Independence of the complex on calcium ion.

Polyacrylamide gels of the complexes were evaluated. Lanes 2-4 show the electrophoretic pattern of covalently linked rTnI-153:TnC complex in the presence of EDTA (a metal chelator). Lane 10 5 shows a non covalently linked rTnI-153:TnC complex in the absence of EDTA. Lanes 6-7 show the electrophoretic pattern of non covalently linked rTnI-153:TnC complex in the presence of EDTA. The disappearance of the upper band of complex is 15 evident in Lanes 6-7, but not in Lanes 2-4. Lane 8 is TnC. See, Figure 7.

Example 5

Preparation of Troponin T Controls.

Example 1 is repeated except that troponin T is substituted for the troponin I. Aliquots of stock solution are diluted to the clinical range of interest from 0.01 ng/mL to 10 ng/mL. The aliquots may be analyzed for immunological activity using a troponin T assay such as that assay commercially available from Boehringer Mannheim.

Example 6

Preparation of a ternary complex of Troponin I, Troponin C and Troponin T.

Example 1 is repeated except that both troponin I and troponin T are added to the activated troponin C. The ternary complex may be observed using electrophoresis run in the presence of EDTA. The troponin I immunological activity may be evaluated as described in Example 1 and the troponin T immunological activity may be evaluated as described in Example 5.

We claim:

1. A composition for use in an assay for the
determination of the presence or concentration of
cardiac troponin I the composition comprising: a
complex of cardiac troponin I or a fragment
thereof and troponin C or a fragment thereof
wherein the cardiac troponin I is covalently
coupled to the troponin C.

2. The composition of claim 1 wherein the
composition is used as a diagnostic calibrator,
control or reference material for cardiac troponin
I.

3. The composition of claim 1 wherein the
composition is immobilized on a substrate wherein
the substrate is used in a competitive assay for
cardiac troponin I.

4. The composition of claim 1 wherein the cardiac
troponin I and troponin C are covalently coupled
using a water soluble carbodiimide.

5. The composition of claim 1 wherein the cardiac
troponin I is recombinant.

6. The composition of claim 5 wherein the cardiac
troponin I is a fragment having substantially the

35

amino acid sequence of a CNBr fragment designated rTnI-153.

7. The composition of claim 5 wherein the cardiac troponin I is a fragment of the general sequence X-A-B-Y wherein X comprises any of amino acids 1-27 of full length cardiac troponin I, A comprises residues 28-69 of full length cardiac troponin I, B comprises amino acid residues 70-90 of full length cardiac troponin I, and Y comprises any sequential amino acid sequence of amino acid residues 91-171 of full length cardiac troponin I.

8. The composition of claim 7 wherein X is amino acid residue 17 of full length cardiac troponin I, and Y is selected from the group consisting of amino acid residues 91-91, 91-101, 91-103, 91-110, 91-117, 91-121, 91-131, 91-140, 91-151, 91-155, 91-161 and 91-171 of full length cardiac troponin I.

9. A composition for use in an assay for the determination of the presence or concentration of cardiac troponin T the composition comprising: a complex of cardiac troponin T or a fragment thereof and troponin C or a fragment thereof wherein the cardiac troponin T is covalently coupled to the troponin C.

10. The composition of claim 9 wherein the composition is used as a diagnostic calibrator, control or reference material for cardiac troponin T.

11. The composition of claim 9 wherein the composition is immobilized on a substrate wherein the substrate is used in a competitive assay for cardiac troponin T.

12. The composition of claim 9 wherein the cardiac troponin T and troponin C are covalently coupled using a water soluble carbodiimide.

13. The composition of claim 9 wherein the cardiac troponin T is recombinant.

14. A composition for use in an assay for the determination of the presence or concentration of cardiac troponin T or cardiac troponin I the composition comprising: a complex of cardiac troponin T or fragment thereof, cardiac troponin I or a fragment thereof and troponin C or a fragment thereof wherein the cardiac troponin T and cardiac troponin I are covalently coupled to the troponin C.

15. The composition of claim 14 wherein the composition is used as a diagnostic calibrator,

37

control or reference material for cardiac troponin I or cardiac troponin T.

16. The composition of claim 14 wherein the composition is immobilized on a substrate wherein the substrate is used in a competitive assay for cardiac troponin I or cardiac troponin T.

17. The composition of claim 14 wherein the cardiac troponin I, cardiac troponin T and troponin C are covalently coupled using a water soluble carbodiimide.

18. The composition of claim 14 wherein the cardiac troponin I is a fragment having substantially the amino acid sequence of a CNBr fragment designated rTnI-153.

19. The composition of claim 18 wherein the cardiac troponin I is a fragment of the general sequence X-A-B-Y wherein X comprises any of amino acids 1-27 of full length cardiac troponin I, A comprises residues 28-69 of full length cardiac troponin I, B comprises amino acid residues 70-90 of full length cardiac troponin I, and Y comprises any sequential amino acid sequence of amino acid residues 91-170 of full length cardiac troponin I.

20. A composition for use in an assay for the determination of the presence or concentration of

38

cardiac troponin T or cardiac troponin I the
composition comprising: a complex of cardiac
troponin T or fragment thereof and cardiac
troponin I or a fragment thereof wherein the
5 cardiac troponin T and cardiac troponin I are
covalently coupled.

1/7

Amino Acid Sequences and comparison of SEQ ID NO. 1 and SEQ ID NO. 2

```

-7          -1 C 1
(r)  M-A-S-M-T-L-W-M-A-D-G-S-S-D-A-A-R-E-P-R-P-A-P-A-P- SEQ ID NO. 1
(h)          A-D-G-S-S-D-A-A-R-E-P-R-P-A-P-A-P- SEQ ID NO. 2

I-R-R-R-S-S-N-Y-R-A-Y-A-T-E-P-H-A-K-K-K-S-K-I-S-A-S-R-K-
I-R-R-R-S-S-N-Y-R-A-Y-A-T-E-P-H-A-K-K-K-S-K-I-S-A-S-R-K-

L-Q-L-K-T-L-L-L-Q-I-A-K-Q-E-L-E-R-E-A-E-E-R-R-G-E-K-G-R-
L-Q-L-K-T-L-L-L-Q-I-A-K-Q-E-L-E-R-E-A-E-E-R-R-G-E-K-G-P-

A-L-S-T-R-C-Q-P-L-E-L-T-G-L-G-F-A-E-L-Q-D-L-C-R-Q-L-H-
A-L-S-T-R-C-Q-P-L-E-L-T-G-L-G-F-A-E-L-Q-D-L-C-R-Q-L-H-

A-R-V-D-N-V-D-E-E-R-Y-D-I-E-A-K-V-T-K-N-I-T-E-I-A-D-L-
A-R-V-D-N-V-D-E-E-P-Y-D-I-E-A-K-V-T-K-N-I-T-E-I-A-D-L-

T-Q-K-I-F-D-L-P-G-K-F-K-R-P-T-L-R-P-V-P-I-S-A-D-A-M-M-
T-Q-K-I-F-D-L-P-G-K-F-K-R-P-T-L-P-P-V-P-I-S-A-D-A-M-M-

Q-A-L-L-P-A-P-A-K-E-S-L-D-L-P-A-H-L-P-Q-V-K-K-E-D-T-E-
Q-A-L-L-G-A-P-A-K-E-S-L-D-L-P-A-H-L-K-Q-V-K-K-E-D-T-E-

K-E-N-R-E-V-G-D-W-R-K-N-I-D-A-L-S-G-M-E-G-R-K-K-K-F-
K-E-N-R-E-V-G-D-W-R-K-N-I-D-A-L-S-G-M-E-G-R-K-K-K-F-
210
|
E-S-P-M-V-H-H-H-H-H-H
E-S
|
209

```

FIG. 1

2/7

SEQ ID:3

A-D-G-S-S-D-A-A-R-E-P-R-P-A-P-A-P-I-R-R-R-S-S-N-Y-R-A-Y-A-
T-E-P-H-A-K-K-K-S-N-I-S-A-S-R-K-L-Q-L-K-T-L-L-L-Q-I-A-K-Q-
E-L-E-R-E-A-E-E-R-R-G-E-K-G-R-A-L-S-T-R-C*-Q-P-L-E-L-T-G-L-
G-F-A-E-L-Q-D-L-C*-R-Q-L-H-A-R-V-D-K-V-D-E-E-R-Y-D-I-E-A-K-
V-T-K-N-I-T-E-I-A-D-L-T-Q-K-I-F-D-L-R-G-K-F-K-R-P-T-L-R-R-
V-R-I-S-A-D-A-M


c*=(cam)=S-carboxyamidomethylcysteine

FIG. 2

SUBSTITUTE SHEET (RULE 26)


3/7

Stability of TnI-153:TnC Complex

**FIG. 3**


4/7

Stability of TnI-153:TnC Complex

**FIG. 4**


5/7

Stability of TnI-153:TnC Complex at 45°C

**FIG. 5**

6/7


Stability of the crosslinked rTnI:TnC Complex at 45°C

**FIG. 6**



7/7

Electrophoretic Gels



Lane: 2 3 4 5 6 7 8

FIG. 7

INTERNATIONAL SEARCH REPORT

International Application No.

PCT/US 98/10513

A. CLASSIFICATION OF SUBJECT MATTER

IPC 6 C07K14/47 G01N33/63 G01N33/53

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC 6 C07K G01N

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
Y	T. KOBAYASHI ET AL: "Extensive interactions between troponins C and I. zero-length cross-linking of troponin I and acetylated troponin C" BIOCHEMISTRY, vol. 34, 1995, pages 10946-10952, XP002078688 see the whole document ---	1-20
Y	J. LESZYK ET AL.: "Cross-linking of rabbit skeletal muscle troponin subunits..." BIOCHEMISTRY, vol. 27, 1988, pages 6983-6987, XP002078689 see the whole document --- -/--	1-20

☒ Further documents are listed in the continuation of box C.

☒ Patent family members are listed in annex.

* Special categories of cited documents

- "A" document defining the general state of the art which is not considered to be of particular relevance
- "E" earlier document but published on or after the international filing date
- "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- "O" document referring to an oral disclosure, use, exhibition or other means
- "P" document published prior to the international filing date but later than the priority date claimed

- "T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
- "X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.
- "3" document member of the same patent family

Date of the actual completion of the international search

25 September 1998

Date of mailing of the international search report

08/10/1998

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentlaan 2
NL - 2280 HV Rijswijk
Tel. (+31-70) 340-2040, Tx. 31 551 epo nl.
Fax: (+31-70) 340-3016

Authorized officer

Cervigni, S

INTERNATIONAL SEARCH REPORT

International Application No.

PCT/US 98/10518

C.(Continuation) DOCUMENTS CONSIDERED TO BE RELEVANT		
Category	Citation of document, with indication where appropriate, of the relevant passages	Relevant to claim No.
Y	EP 0 743 522 A (PASTEUR SANOFI DIAGNOSTICS) 20 November 1996 see the whole document ---	1-20
A	WO 94 27156 A (FORTRON BIOSCIENCE INC) 24 November 1994 ---	
A	MORJANA N A ET AL: "BIOCHEMICAL AND IMMUNOLOGICAL PROPERTIES OF A CYANOGEN BROMIDE FRAGMENT OF HUMAN CARDIAC TROPONIN I" FASEB JOURNAL, vol. 10, no. 6, 30 April 1996, page A-1295 XP002028680 ---	
A	C.S. FARAH ET AL.: "The troponin complex and the regulation of muscle contraction" FASEB JOURNAL, XP002078690 see the whole document -----	

INTERNATIONAL SEARCH REPORT

Information on patent family members

International Application No

PCT/US 98/10518

Patent document cited in search report		Publication date	Patent family member(s)	Publication date
EP 0743522	A	20-11-1996	FR 2734267 A	22-11-1996
			AT 169740 T	15-08-1998
			CA 2176655 A	17-11-1996
			DE 69600512 D	17-09-1998
			JP 9021804 A	21-01-1997
<hr/>				
WO 9427156	A	24-11-1994	AU 689777 B	09-04-1998
			AU 6834894 A	12-12-1994
			CA 2162121 A	24-11-1994
			EP 0699306 A	06-03-1996
			JP 9503050 T	25-03-1997
			US 5756682 A	26-05-1998
<hr/>				

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.

This Page Blank (uspto)